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GB User Instructions

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boso
medicus
family

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Package

- 1 Blood pressure unit
boso-medicus family
- 1 Cuff CA01
- 1 Soft case
- 2 Blood pressure record cards
- 1 Certificate of warranty
- 4 Batteries LR 6
- 1 User instructions

Preliminary Notes

Your **boso-medicus family** unit is a state of the art blood pressure computer. It works on the **oscillometric measuring principle**. The oscillations created by the pulse and transmitted via the cuff are stored and evaluated by a micro processor. Therefore there is no need for a microphone, the positioning of which is vital if one is to get reliable readings.

The boso-medicus family conforms to the requirements of the European Directives based on the Law governing Medical Products (CE mark), and to the European Standard EN 1060, Part 1, on „Non-Invasive Blood Pressure Instruments – General Requirements“ and Part 3 „Additional Requirements for Electronic Blood Pressure Systems.“

Regular checks on the accuracy of the measurement system should be carried out in accordance with the valid directives of the individual countries, governing this issue.

The instructions for re-calibration checks are on page 33 of this booklet.

Blood Pressure (BP) Values

To determine correct blood pressure, the following two values have to be ascertained:

- The **Systolic or upper blood pressure**. This is created when the heart muscle contracts and blood is being forced into the blood vessels.
- The **Diastolic or lower blood pressure**. This is created when the heart muscle extends and once again fills with blood.

Measurement values are expressed in **mmHg (mm mercury column)**.

The World Health Organization (WHO) has set the following guideline for the assessing of blood pressure values:

	Systolic	Diastolic
Too high	over 140 mmHg	over 90 mmHg
Normal – borderline	130 to 139 mmHg	85 to 89 mmHg
Normal	120 to 129 mmHg	80 to 84 mmHg
Optimum	up to 119 mmHg	up to 79 mmHg

Whether medical treatment is necessary depends not only on the blood pressure but also on the risk profile of the individual patient. Please contact your GP if one of the values (SYS / DIA) constantly exceeds the limit (> 140 / > 90).

Important Observations to be followed when Self Testing

- Always take your blood pressure in a calm and relaxed state. It is recommended to do two checks daily, – on getting up in the morning and again at night before retiring.
- The blood pressure is best taken sitting down. Sit comfortably in an armchair or at a table where you can easily support your left arm, holding it slightly bent.
- The measurement is taken on the free left upper arm. Please ensure that there is no restriction from a tight, pushed-up sleeve. It may be necessary to remove the garment.
- During measurement the arm has to be kept absolutely still.
- Between successive measurements an interval of at least two minutes must be allowed.
- Differences in blood pressure are quite normal. Even in closely repeated measurements there can be marked differences. One single measurement, or measurements taken at irregular intervals, do not give any clear indication of the real blood pressure. An accurate evaluation is only possible when measurements are taken on a regular basis and for the established values to be entered in the blood pressure check card.
- Self testing does not mean self treatment. Never change the doses of any medication as prescribed by your doctor.
- Irregular heart beats can influence the sensitivity of the unit and result in inaccurate measurement values being displayed. Adverse results can also be experienced by patients with a pacemaker who have a weak pulse. The BP unit itself has however no detrimental effect on the pacemaker.
- Using your blood pressure monitor in the vicinity of strong electro-magnetic fields, (ie mobile telephones) can also cause faulty readings.

Starting up

Fitting the Batteries

The battery compartment is located on the underside of the unit. Insert batteries as shown in the illustration.


Please note: If batteries are inserted incorrectly the unit does not function and this can lead to warming and leakage.

- Only use high quality leak proof batteries of the correct specifications (see Section Technical Data) and never mix used with new batteries or batteries of different manufacture.
- Remove batteries if the unit is not to be used for any length of time.
- The unit shows the following function symbols:

 battery fully loaded

 battery partially loaded

 replace battery soon

 flashing! No further measurements possible; replace battery.

- **Care for your environment**

Used batteries and recharging units do not belong into the household dustbin.

They can be deposited at special refuse centres for used batteries. Contact your townhall for information.

- **Change batteries only when the unit is switched off!**

The measurement values will be deleted if the batteries are not replaced within 30 seconds. With the removal of the batteries while the unit is switched on, the memory store is erased immediately.

How to take your own Blood Pressure

Fitting the cuff


- Open cuff into a ring shape and pull onto your upper left arm until the lower cuff edge lies 2 – 3cm above the elbow ensuring that the special marking is positioned directly on the artery. The metal ring must never come to lie on the artery as this may cause a false reading.
- Gently pull the cuff end that feeds through the metal ring and place it up and outside around your arm. Close velcro fastening with a light pressure. The cuff must not be too tight. There should roughly be space for two fingers between arm and cuff.
- Note on cuff sizes:
This unit must only be used with the following cuff types:

Type	Size of Arm Rubber Bag	Circumference	Order Number
CA01	13 x 25 cm	22 – 32 cm	143-4-750*
CA02	14 x 30 cm	33 – 41 cm	143-4-757

* comes as standard

Important :

A measurement should not take longer than 2 minutes. In case of an error the measurement can be interrupted at any time by opening the velcro fastening of the cuff.

Taking your Blood Pressure


- Attach cuff the BP monitor by fitting the plug firmly into (click) the opening at the left side of the unit.


- Place the arm with the cuff on the table holding it relaxed and slightly bent so that the cuff is line with your heart.

- Choose a symbol for your measurement:

- Press the button with your symbol briefly. The symbol button is illuminated in green colour. All digits and function symbols in the display are momentarily lit up to indicate that the unit is now ready for use.

- **Note for guest mode:**

If the measurement shall not be stored (guest mode), the measurement is started by pressing both buttons briefly at once. No button is illuminated.

- At the first measurement after fitting the batteries the unit automatically inflates the cuff to approx. 160 mmHg. The rising cuff pressure is displayed by the moving digits. Should the pressure not be sufficient the unit increases this automatically by approx. 40 mmHg.

- During the next measurement the inflation height orientates itself against the result of the last measurement (**memomatic** function: inflation value = Systole of last measurement + 40 mmHg).
- When the correct pressure is reached the pump switches off and the air is automatically released from the cuff.
- While the blood pressure is being taken the drop in cuff pressure is again digitally displayed and symbol ♥ is lit up. Remember to remain perfectly still and not to speak.

Display of Measured Values

- As soon as the first pulsations are registered, symbol  flashes in Pulse rhythm.
- After completion of a measurement the built-in magnetic valve opens automatically to allow rapid deflation of the cuff.
- For each symbol  and  the unit has an individual memory store for 30 measurements each.
- The measurement is stored automatically (Exception: guest mode, see page 25). As soon as your memory store is full (30 measurements maximum) the oldest measurement gets deleted and the latest one is stored.
- Invalid measurements are not stored in the memory.
- If symbol  appears following a completed measurement, it is recommended that the measurement is repeated while holding the arm motionless.
A renewed display of symbol  following the repeat measurement with the arm being held motionless, can be an indication of irregular heart beat. Discuss this matter with your GP on your next visit.
- The unit switches off automatically after approximately one minute. To repeat a measurement press your symbol button again briefly.
- If no further measurement is required, remove cuff.


Interruption of a measurement :

A measurement can be interrupted at any stage by simply pressing your symbol button again. This triggers the automatic air release from the cuff (For further measurement press your symbol button again).

How to record blood pressure values in your record card and recall the average value (including evaluation according WHO)

- Note the measured values in your blood pressure record card (two cards are included).
- Enter each individual measurement.
- Recall the individual average value after 30 measurements and note it in its designated column of the BP record card. The average value is recalled by pressing the relevant symbol button for about 3 seconds while the unit is switched off. The average value of all stored values (systole and diastole) appears. In the „pulse“ display the number of retained measurements flashes.
- The symbol button of the displayed value is illuminated in green colour.
- If there is no measurement stored in the memory, the display for the number of measurements taken flashes „000“. The systole and diastole displays each show „0“.
- Only the left side of the display appears an index relating to the average value. The evaluation scale immediately shows you where your blood pressure value ranges, according to the WHO classification (see page 20).


Evaluation in accordance with WHO classification

- Too high
- Normal – borderline
- Normal
- Optimum

How to recall stored measurement values

- The latest measurement stored will be indicated immediately after the average value has been displayed.

The number of the measurement appears in the display. After 3 seconds the corresponding measurement value (systole, diastole, pulse) appears automatically.

- The symbol button of the displayed value is illuminated in green colour.
- All stored measurements are displayed in the same manner automatically.
- The unit can be switched off by simply pressing the symbol button.
- If no button is pressed, the unit switches off automatically after 5 seconds.
- **Change batteries only when the unit is switched off!**
The measurements values will be deleted if the batteries are not replaced within 30 seconds. With the removal of the batteries while the unit is switched on, the memory store is erased immediately.

Systole
mmHg

Diastole
mmHg

Puls/min.

Mains Operation


- A mains socket connector is located on the rear side of the unit. For operating the unit from the mains use only the **boso power supply pack (order number 410-7-150)**. This power supply pack produces a rectified output of the correct polarity. Other commercially available power supply units can cause damage to the electronic components resulting in the guarantee becoming void.

Important: If no batteries are in the unit, the memory store is cleared as soon as the power supply lead is disconnected. Insert batteries in order to save the memory store.

Switch off the unit, disconnect the connector from the unit and then the adapter from the wall socket.

Error Messages

- Should there be any problem while taking your blood pressure that impedes correct monitoring, the display will show an error message instead of a BP reading.
- Meaning of the error messages:

Err: cuff pressure varied while blood pressure is taken
→ hold arm still

– or –

no valid pulse readings could be obtained
→ check position of cuff

– or –

Systolic – Diastolic pressure ≤ 10 mm/g
→ check position of cuff

Err Incorrect inflation

CUF: → cuff possibly not firmly enough positioned

«♥» : Irregular puls or incorrect taking of blood pressure (i.e. movement during measurement)
→ Repeat measurement holding the arm absolutely still. If the «♥» symbol is displayed again, it can be that an irregular heart beat is present. Consult your GP.

Cleaning of BP Unit and Cuff

- To clean your unit use a soft dry cloth only.
- Small stains can be removed with a proprietary stain remover.

Guarantee information

For details on the conditions of guarantee, please refer to the certificate of warranty supplied.

A claim for guarantee can only be allowed if the Certificate of Warranty, completed and stamped by the dealer, is enclosed with the instrument.

For both guarantee and repair work, please send the instrument carefully packed and postage paid either to your authorized dealer or directly to:

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Technical Data:

Measurement system: Oscillometric

Measurement range: 40 – 240 mmHg, 40 – 200 Puls/min.

Cuff pressure: 0 – 320 mmHg

Memory store: 2 x 30 measurements

Display: LCD

Working conditions:

Environmental temperature: Room temperature 10 – 40°C

Relative humidity 10 – 85 %

Storing conditions:

Environmental temperature -5°C to +50°C

Relative humidity 85 % maximum

Power supply:

DC 6 V (batteries 4 x 1,5 Mignon IEC LR 6, alkaline manganese)

Alternative special option: Main unit DC 6 V, Order-No. 410-7-150 (Polarity: outside MINUS, inside PLUS) 

Battery check: Symbol display in LCD window

Weight: 300 g without batteries

Dimensions (B x H x T):

150 mm x 65 mm x 115 mm

Classification: Class II ()

Typ BF ()

Clinical Test (DIN 58130):

Accuracy complies with EN 1060 part 3

Maximum deviation of cuff pressure: ± 3 mmHg

Maximum deviation of pulse rate: ± 5 %

The right to alter technical specifications is reserved.

A) Function testing

A function test can only be carried out on a person or with a suitable simulator.

B) Testing of the Density of Pressure Circuit and Divergence in the Pressure Indication

Please note:

- a) When in the measurement mode pressure is increased to above 320 mmHg, the quick release valve is triggered and the pressure circuit opened. When during testing pressure is increased to above 320 mmHg the measurement display will flash.
- b) Prior to testing the connector piece has to be removed from the cuff tube. Then the short connecting head of the connector plug – up to now attached to the unit and its corresponding socket at the unit – has to be attached to the cuff tube, and the long connecting head of the connector piece to be inserted into the unit.
Lastly the pump ball is to be fitted into the pressure circuit.

Testing

- 1.) Remove batteries.
- 2.) Hold symbol button  pressed down and insert batteries.
- 3.) Release symbol button .
The unit is now ready for testing. The displays of SYS, DIA and PUL now show the actual pressure. The symbol button  is illuminated in green colour.
- 4.) Testing for divergence of Pressure indication and Density of Pressure circuit (setting time for cuff at least 30 seconds) can now be carried out in the usual manner.
- 5.) Switch off unit by pressing the symbol button  again and re-establish all normal connections.

C) Safety Seal

As a measure to safeguard the correct calibration of the BP unit the upper and lower parts of the housing can be connected with a safety seal. Alternatively, a hole drilled for the connecting screw at the back of the unit can be sealed for safety.

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**BOSCH + SOHN GmbH u. Co. KG · Fabrik mediz. Apparate