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Ⓒ **Instructions for use**

Ⓕ **Mode d'emploi**

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boso

medilife

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For details on the conditions of guarantee, please refer to the guarantee card supplied.

A claim for guarantee can only be allowed if the Certificate of Warranty, completed and stamped by the dealer, is enclosed with the instrument.

For both guarantee and repair work, please send the instrument carefully packed and postage paid either to your authorized dealer or directly to:

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Important notes	55	Recalling stored readings	64
Description of unit	56	Error indicator	65
Principle of operation	56	Cleaning and care of the unit	66
Information about blood pressure	57	Environmental safety notes	66
Blood pressure readings	58	Technical data	67
Preliminary remarks on measuring your own blood pressure	59	Technical testing procedures	68
Use and function	60		
Inserting the batteries	60		
Preparing to measure blood pressure	61		
How to wear the unit	61		
Measuring your blood pressure	62		
A step-by-step guide	62		
memomatic function	62		
Reading the display	63		

- Do not expose the unit to extreme temperatures and avoid exposure to prolonged direct sunlight. Guard against mechanical shocks (do not let the unit fall or be struck) and do not allow the unit to become soiled or wet.
- To prevent possible damage to the unit from leaking batteries, remove the batteries if the unit is not going to be used for a longer period.
- Operating the unit within the vicinity of strong electro-magnetic fields (e.g. radiometric applicators, mobile phones) may cause the unit to malfunction.
- The cuff of the sphygmomanometer is suitable for a wrist circumference between 13.5 and 21.5 cm.
- Take care to ensure that the cuff is not damaged by sharp or pointed objects (needles, scissors, etc.)

Principle of operation

Your boson sphygmomanometer has been designed and constructed on the basis of the latest technological developments. It operates on the oscillometric principle. The pressure variations (oscillations), resulting from the pulse waves and transmitted by the cuff, are stored and evaluated by a microprocessor.

The major advantage of this method is that no microphone is required. Normally, the reliability of a blood pressure measurement is strongly dependent on the accuracy with which the microphone is positioned when the cuff is put into place.

This sphygmomanometer complies with the European directives underlying the guidelines for medical products (CE symbol), as well as the European Standard EN 1060 Part 1: 'Non-invasive Sphygmomanometers - General Requirements' and Part 3: 'Supplementary Requirements governing electronic blood pressure measurement systems'. Regular calibration testing is to be carried out in accordance with the regulations existing in individual countries.

The technical testing procedures are detailed on page 68.

The heart, which can be thought of as a continually active 'pump' driving the blood under pressure through our body's circulatory system, ensures that our organs receive an adequate supply of blood vital to our survival.

Blood pressure depends on a number of factors, including the strength of the heart beat and the elasticity and diameter of the blood vessels. The blood pressure which arises when the heart muscle contracts (systolic pressure) is higher than that present when the heart muscle relaxes (diastolic pressure).

Blood pressure does not remain constant but is continually changing. In the morning the blood pressure is lower than in the evening and it is at its lowest during sleep. It is also influenced by eating, smoking, stress and bodily exertion.

Transient increases in blood pressure are no cause for concern. However, if the blood pressure is permanently above the normal value, one speaks of high blood pressure (hypertension).

High blood pressure is not normally noticeable immediately by the occurrence of particular complaints. If, though, it goes unnoticed or untreated, it can lead to serious health problems. For this reason, regular blood pressure measurements as well as diagnosis and treatment by a doctor are so important.

A reliable assessment of blood pressure is only possible if you measure blood pressure regularly under comparable conditions and record the readings in your blood pressure log book.

Measuring your own blood pressure is not, however, a therapy. Never alter on your own accord any medication dosage that your doctor has prescribed.

Blood pressure readings

To calculate the blood pressure, two values must be measured:

1. The systolic (upper) blood pressure value.

This arises when the heart contracts and blood is forced into the blood vessels.

2. The diastolic (lower) blood pressure value.

This is the value present when the heart muscle expands and fills again with blood.

Blood pressure readings are expressed in mmHg (millimeters of mercury).

The World Health Organization (WHO) has established the following parameters for the evaluation of blood pressure values.

Normal blood pressure:

Systolic – up to 140 mmHg

Diastolic – up to 90 mmHg

Please consult your doctor if your own results consistently exceed even only one of the above figures.


- In certain cases, a different blood pressure reading may be recorded at the wrist to that measured on the upper arm.

We therefore recommend that comparison measurements are occasionally made on the upper arm.

- Pregnant women, diabetics and persons suffering from liver disorders, cardiac dysrhythmia (heart rhythm disturbance) or those with a cardiac pacemaker should consult their doctor before measuring their own blood pressure, as in such cases, altered or falsified blood pressure readings may occur. The sphygmomanometer does not have any influence on the pacemaker.
- The measurement should always be performed when relaxed and when sitting. Measurement should be carried out on the bare left wrist.
- The arm must not be constricted as a result of clothing that has been pushed back or rolled up, and the arm should be kept still throughout the measurement.

- **When measuring, the wrist must be held at the same height as the heart!**

- It is recommended that the blood pressure is measured twice a day. In the morning after getting up and in the evening when the body is relaxed.
- Measurement on the wrist is not suitable for persons suffering from circulatory disturbances or with severe arteriosclerosis (thickening of the artery walls). In these cases, blood pressure measurement should be performed on the upper arm.
- Variations in blood pressure are perfectly normal. Even a repeat measurement can produce noticeable differences. One-off or irregular measurements are therefore of little value when trying to assess the actual blood pressure. An interval of at least two minutes should be observed, before repeating a measurement.

Inserting the batteries

The battery compartment is located on the upper side of the casing.


To open, push the catch over to the left.

Insert two 1.5 volt micro-cells (LR03) in the battery compartment as illustrated on the sleeve.

Only use the specified batteries. Never mix old and new batteries or different makes. Then close the battery compartment again.

Empty batteries must not be disposed of as household waste (see notes on environmental notes, page 66).

The batteries require replacing when the symbol  is illuminated in the display.


How to wear the unit

Wrap the cuff of the unit around your left wrist. Ensure that the display is on the underside of your arm (see drawing).


Close the cuff with the velcro fastener. The cuff should fit snugly around the wrist.

To ensure that the results are as accurate and clear as possible, sit down and relax for a few minutes before starting the measurement.

Hold your arm in a relaxed position to your body (see drawing). In this position, your wrist joint with the measuring unit is on a level with your heart.

Note on cuff size: The unit must only be used in combination with the CW 31 cuff (wrist circumference 13.5 to 21.5 cm).

A step-by-step guide

- Press the START button. All segments of the numerical display and all function symbols will be illuminated for a short period.
- At the first measurement after insertion of the batteries the unit automatically inflates to ca. 180 mmHg. The increasing sleeve pressure is


shown digitally.

If the reached pressure level is too low, the unit automatically pumps up by approx. 50 mmHg. Successive measurements adjust the inflation level according to the result of the previous measurement – **memomatic** function: Inflation level systolic value of the foregone measurement + 50 mmHg.

- Once inflated, the pump is switched off and the air in the cuff escapes automatically.
- If not already the case, you should now keep your arm completely still and not speak.

- As soon as the unit registers the first beats, the ♥ symbol blinks in time with the pulse rhythm and the falling pressure in the cuff is displayed
- After the measurement, the built-in solenoid valve opens automatically to rapidly release the remaining pressure in the cuff.

The measurement readings are shown in the display:


- The unit switches off automatically after approximately one minute. To repeat a measurement, press the START button again. For **boso-medilife PC** the unit changes to the time display mode.

Aborting a measurement

- A measurement can be aborted at any time by pressing the START button. The cuff is then automatically deflated. In case of error, you can always remove the unit at any time by simply opening the cuff.

Recalling stored readings

- The unit stores the 30 most recent measurements.
- Faulty measurements are not saved.


In front of the first saved measured value, the display shows $\square 01$.
In front of the second saved measured value, the display shows $\square 02$ etc.

- All measurements in the memory are now shown one after the other for approx. 3 seconds. (If there is no measurement in the memory, then the display shows a „0“ in each case for systole, diastole and pulse).
- Press the START key again to interrupt the display of the saved values and start a new measurement.

- Up to 30 measurements are saved. Once the measured value memory is full, the oldest measurement is then deleted and the latest measurement saved.

To access saved measurements, when you switch the unit on, hold the START key until the last saved measurement appears in the display (symbol **M**).

Error indicator

If a problem should arise during a measurement which prevents a usable measurement from being made, an error number is displayed in place of the usual measurement result.

Error number key:

Err: Inflation pressure in-adequate after three attempts to inflate.
– Keep arm still.

No pulses were detected that could be evaluated.
– Check the positioning of the cuff.

Err CUF: Error during inflation.
– Cuff is possibly not tight enough.

Err B,
Err C,
Err M,
Err S: Please send the unit to the manufacturer's to be checked (adress see page 53).

Cleaning and care of the unit

Please use only a soft, dry cloth to clean the unit.

Never use thinners, petroleum spirit or other aggressive cleaning agents.

Environmental safety notes

Packaging materials or a unit that is no longer of use should be recycled and not simply discarded.

Empty batteries and storage batteries should not be thrown away with household waste. They can be handed in at a battery collection point or disposed of as special waste.

Please contact your local authority for details.

Technical data

Measurement principle: oscillometric

Measurement range:

40 to 260 mmHg

40 to 200 pulses/min

Display: LCD

Operating conditions:

ambient temperature: 10 °C – 40 °C

relative humidity: 10 – 85 %

Transport and storage conditions:

Temperature -5 °C – 50 °C

relative humidity: max. 85 %

Power supply:

DC 3V (2 x 1.5 V microcells alkali manganese, IEC LR 03)

Battery check: indicator symbol in display

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diastolic blood pressure: ±5.53 mmHg.

Accuracy:

maximum measuring deviation of the sleeve pressure: ± 3 mmHg

Error margin of the pulse frequency measuring device is ± 5 %

Technical specifications subject to change without prior notification.

A) Functional testing

Functional testing of the unit can only be performed on human subjects or using a suitable simulator.

B) Leak testing of the pressure circuit and error testing of the display**Note:**

- a) For testing to be carried out, the cuff must be removed from the unit. On the side of the sleeve (at the top and bottom of the unit) there are two clips. To remove the sleeve, push the clips out with a pointed tool (e.g. tweezers). Once the sleeve is removed, the air connection is visible.
- b) To perform leak testing of the pressure circuit, a tee connector must be used to connect the rolled-up cuff and a pump bulb to the pressure circuit.

- c) When the pressure in the measuring mode exceeds 320 mmHg, the quick drain valve opens. When the pressure in the test mode exceeds 320 mmHg, the display flashes.

Testing

1. Connect the unit to be tested with the pressure generator (e.g. pump bulb).
2. Press the START button briefly. The unit inflates the pressure circuit to about 180 mmHg.

3. Immediately after inflating, raise the pressure rapidly to around 320 mmHg so that the quick-action air release valve triggers opening the pressure circuit.
4. Wait until the pressure has fallen to 0 mmHg and then raise the pressure to about 100 mmHg. The unit is now in test mode, the current pressure values are then shown in the display regions: SYS, DIA and PUL.
5. Carry out error testing of the pressure reading in the normal manner.

6. Perform leak testing in the usual fashion bearing in mind the settling time of two minutes.
7. Switch off unit and reconnect the cuff to the unit.

C) Security

For security purposes the two parts of the housing are joined together with a security label.

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**BOSCH + SOHN GmbH u. Co. · Fabrik mediz. Apparate